

Data Sheet



QT 70 SDA

70% Specially Denatured Ethanol with IPA

QT 70 SDA is a general-purpose disinfectant for use in controlled environments. QT 70 SDA is a blend containing 70% Denatured Ethanol with IPA and 30% USP WFI filtered to 0.22 micron and Gamma Irradiated to SAL of 10^{-6} .

Each bottle is lot number and expiry date coded for complete traceability and straightforward record keeping. Bottles come fully assembled with a non-venting, ergonomic trigger sprayer with mist. The closed system (Bag in a Bottle) spray assembly maintains sterility throughout the products use. All sterile product is individually double-bagged in heat sealed poly bags and packaged in cases with individually sealed double case liners.

Features

- Contains 70% Specially Denatured Alcohol and 30% WFI
- 0.22 micron filtered
- Gamma irradiated to SAL of 10^{-6} at 25-45 kGy
- Double-bagged and double sealed polybag case liner
- Lot specific records for each shipment:
 - Certificate of Irradiation
 - Certificate of Analysis
 - Indelible lot & expiry markings on each container

Benefits

- Removes most surface contaminants
- Low residue upon evaporating/no rinsing
- Pre-attached and ready to use trigger sprayer

Applications

- Surface cleaning
- Equipment cleaning
- Gloved hand cleaning in sterile areas
- Wipe down for pass through to controlled environments



* Non-Sterile product will **NOT** include Gamma Indicators or Double-bagged.

Part No.	Description	Quantity	Dimension	Weight	HS Code	UN Code	Shelf Life
QT2512-16	70% SDA, 30% WFI Bag in Bottle – 473mL – Sterile	12 x 500mL/case	350mm (L) x 280mm (W) x 248mm (H)	6.5 kg	3808 94 9000	UN1170	36 months
QT2512-32	70% SDA, 30% WFI Bag in Bottle – 946mL – Sterile	6 x 1L/ case	305mm (L) x 215mm (W) x 320mm (H)	6.5 kg	3808 94 9000	UN1170	36 months
QT2513-16	70% SDA, 30% WFI Trigger Spray – 473mL - Sterile	12 x 500mL/case	350mm (L) x 280mm (W) x 248mm (H)	6.5 kg	3808 94 9000	UN1170	36 months
QT2513-16NS	70% SDA, 30% WFI Trigger Spray – 473mL	12 x 500mL/case	350mm (L) x 280mm (W) x 248mm (H)	6.5 kg	3808 94 9000	UN1170	36 months
QT2513-32	70% SDA, 30% WFI Trigger Spray – 946mL- Sterile	6 x 1L/ case	305mm (L) x 215mm (W) x 320mm (H)	6.5 kg	3808 94 9000	UN1170	36 months
QT2513-32NS	70% SDA, 30% WFI Trigger Spray – 946mL	6 x 1L/ case	305mm (L) x 215mm (W) x 320mm (H)	6.5 kg	3808 94 9000	UN1170	36 months
QT2514-128	70% SDA, 30% WFI Screw Cap – 1 gal – Sterile	4 x 1 gal/case	410mm (L) x 260mm (W) x 300mm (H)	15.0 kg	3808 94 9000	UN1170	36 months
QT2514-128NS	70% SDA, 30% WFI Screw Cap – 1 gal	4 x 1 gal/case	410mm (L) x 260mm (W) x 300mm (H)	15.0 kg	3808 94 9000	UN1170	36 months
QT2530-9x20	Presat Polyester Wipes 70% SDA, 30% WFI, 9”x 9” – Sterile	10 pouches/case	425mm (L) x 279mm (W) x 318mm (H)	8.0 kg	3808 94 9000	UN3175	24 months
QT2530-9X20NS	Presat Polyester Wipes 70% SDA, 30% WFI, 9”x 9”	10 pouches/case	425mm (L) x 279mm (W) x 318mm (H)	8.0 kg	3808 94 9000	UN3175	24 months
QT2530-12x20	Presat Polyester Wipes 70% SDA, 30% WFI, 12”x12” – Sterile	10 pouches/case	425mm (L) x 279mm (W) x 318mm (H)	8.0 kg	3808 94 9000	UN3175	24 months
QT2530-12X20NS	Presat Polyester Wipes 70% SDA, 30% WFI, 12”x12”	10 pouches/case	425mm (L) x 279mm (W) x 318mm (H)	8.0 kg	3808 94 9000	UN3175	24 months
QT2531-911x20	Presat Non-Woven Cellulose Polyester Wipes 70% SDA with IPA, 30% WFI, 9”x11” – Sterile	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months
QT2531-911x20NS	Presat Non-Woven Cellulose Polyester Wipes 70% SDA with IPA, 30% WFI, 9”x11”	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months
QT2532-911x20	Presat Non-Woven MBPP Wipes 70% SDA with IPA, 30% WFI, 9”x11” – Sterile	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months
QT2532-911x20NS	Presat Non-Woven MBPP Wipes 70% SDA with IPA, 30% WFI, 9”x11”	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months
QT2515-11	70% SDA, 30% WFI BOV Aerosol – Sterile	11oz x 24 cans/case	420mm (L) x 300mm (W) x 245mm (H)	10.5 kg	3808 94 2000	UN1950	36 months
QT2515-11NS	70% SDA, 30% WFI BOV Aerosol	11oz x 24 cans/case	420mm (L) x 300mm (W) x 245mm (H)	10.5 kg	3808 94 2000	UN1950	36 months

PRODUCTS AVAILABLE STERILE OR NON-STERILE.

PMA QuantumTec Life Science Solutions

- No.11, Lintang Beringin 3, Diamond Valley, 11900 Penang, Malaysia.
- Website: www.quantumtec.com.my
- Tel: 604-626 5518
- Email: pmasales@pma-asia.com