



Data Sheet

QT PRE-SATURATED NON-WOVEN CP WIPES

QT's Hydroentangled Cellulose/Polyester wipers are a cost-effective cleanroom wiper providing low levels of particles and fibers at an economical price. These highly absorbent wipers offer good wet strength and can be used for many general cleaning applications. The perfect wiper for sanitizing articles prior to pass through, routine cleaning and wiping down of tools, instruments and other equipment. QT's Pre-Saturated Hydroentangled Cellulose/Polyester wipers are available with 70% IPA, 70% SDA and **RTU Detergent.**

Non-Woven CP Wipes

The Nonwoven cellulose polyester is a hydroentangled 55% cellulose and 45% polyester nonwoven wiper designed for applications that require the unique combination of superior strength and exceptional absorbency. This wiper is a perfect wiping material for cleaning and general-purpose wiping.

Features

- Made of hydroentangled 55% cellulose and 45% polyester
- Cleanroom processed and packaged to reduce particles
- Combines the strength of polyester with the absorbency of long cellulose fibres
- Smooth, highly absorbent and durable
- Works well with IPA and other solvents
- 9"x 11"

Applications

- Recommended for use in Class 100 100,000 cleanroom environments ٠
- General purpose wiping of surface areas, including equipments
- Excellent for use in Pharmaceutical, Biotech & medical device manufacturing markets
- Acceptable for use in aseptic areas

Technical Data

Basis Weight	68 ±4 gsm	
Thickness	0.27 - 0.34 mm	
Absorbency Capacity (IEST-RP-CC004.3, Sec 8.1)	≥ 250 mL/m²	
Absorbency Rate (IEST-RP-CC004.3, Sec 8.2)	< 1.0 sec	
Non-volatile Residue, NVR (IEST-RP-CC004.3, Sec 7.1.2)		
In WFI water In IPA	0.021 g/m² 0.007 g/m²	
Organics by FTIR (IEST-RP-CC004.3, Sec 7.2.1)	No Silicone oil, Amide, DOP detected	
Particles by LPC (IEST-RP-CC004.3, Sec 6.1.4)	6.1.4) 75.0 x 10 ⁶ counts/m ²	
Fibre, Particles >100µm (IEST-RP-CC004.3, Sec 6.1.4, 6.2.2.1, 6.2.2.2)	86,181 counts/pc	
Particles by APC sizes@≥0.5μm (IEST-RP-CC003.2, Helmke drum @ 10rpm)		
lonics, IC (IEST-RP-CC004.3, Sec 7.2.2)		Quantu
Chloride (Cl)	68.79 ppm	No. 12, 1971 No. 13, 1971 No. 13, 1974 No. 10, 10, 10, 10, 10, 10, 10, 10, 10, 10,
Sodium (Na)	33.89 ppm	A MERCEN AND
Potassium (K)	7.051 ppm	And and a second a
Magnesium (Mg)	4.63 ppm	ALL DE LE DE
Calcium (Ca)	17.37 ppm	





Note:

Data shown are typical values and should not be used as specifications.

PMA QuantumTec Life Science Solutions

- No.11, Lintang Beringin 3, Diamond Valley, 11900 Penang, Malaysia. Tel: 604-626 5518
- Website: www.quantumtec.com.my

- Email: pmasales@pma-asia.com







QuantumTec offers a complete line of Pre-Saturated wipers in many sizes, packaging quantities and substrates.

All QT Pre-Saturated Wipers are available in both Sterile & Non-Sterile variants.

Each Peel & Reseal pouch is lot numbered and expiry date coded for complete traceability and straightforward record keeping.

All Sterile product is individually double-bagged in heat sealed poly bags and packaged in cases with individually sealed double case Liners.

Sterile wipers are Gamma irradiated to SAL of 10⁻⁶ at 25-45 kGy.

Product Specific Technical Data Sheets are available for Pre-Saturated wipers listed below.

Part No.	Description	Quantity	Dimension	Weight	HS Code	UN Code	Shelf Life
QT1531-911x20	Presat Non-Woven Cellulose Polyester Wipes 70% IPA, 30% WFI, 9"x11" – Sterile	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	2905.12.0000	UN3175	24 months
QT1531-911x20NS	Presat Non-Woven Cellulose Polyester Wipes 70% IPA, 30% WFI, 9"x11"	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	2905.12.0000	UN3175	24 months
QT2533-911x20	Presat Non-Woven Cellulose Polyester Wipes 70% SDA with Methanol, 30% WFI, 9"x11" – Sterile	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months
QT2533-911x20NS	Presat Non-Woven Cellulose Polyester Wipes 70% SDA with Methanol, 30% WFI, 9"x11"	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months
QT2531-911x20	Presat Non-Woven Cellulose Polyester Wipes 70% SDA with IPA, 30% WFI, 9"x11" – Sterile	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months
QT2531-911x20NS	Presat Non-Woven Cellulose Polyester Wipes 70% SDA with IPA, 30% WFI, 9"x11"	10 pouches/case	290mm (L) x 190mm (W) x 250mm (H)	3.70 kg	3808 94 9000	UN3175	24 months

PRODUCTS AVAILABLE STERILE OR NON-STERILE.

PMA QuantumTec Life Science Solutions

- No.11, Lintang Beringin 3, Diamond Valley, 11900 Penang, Malaysia. Tel: 604-626 5518
- Website: <u>www.quantumtec.com.my</u>

- Email: pmasales@pma-asia.com